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Join the Care Connection money saving offers, products to win free products and more with your symptoms.

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OPERATION/MOVEMENT

- 1) Place your hands on the handle bars.
- 2) Place the knee of your injured foot on the pad so that the top of your boot or cast is in the middle of the platform. Your lower leg is then resting on the pad.
- 3) Propel yourself forward or backward with your pushing foot, keeping the knee of your injured foot centered on the pad.
- 4) Move slowly so that you can stop safely using your pushing foot.
- 5) The hand brake may be used if needed but is primarily designed to maintain a stopped stance. The braking wheel must be in contact with the floor, and the brake works best with weight on the knee walker.

DO NOT ATTEMPT TO ADJUST OR OPERATE KNEE WALKER WITHOUT READING ALL INSTRUCTIONS CAREFULLY.

Ask your health professional or knee walker distributor to provide instruction and a demonstration if needed. Failure to properly engage the quick-release clamp may cause the steering column to move while the knee walker is in use and may cause you to lose control. Be sure to check that the cotter pins, clevis pins, nuts, bolts and quick-release clamps are secure and tight before each use or injury may result. Do not use the knee walker if any pins are missing. Use of the hand brake while in motion may cause an abrupt stop.

This is a walking aid only. Do not turn while coasting or if one foot is not on the ground. Do not use the knee walker on or near stairs or while using mind-altering drugs or alcohol or while suffering from dizziness. Do not use the knee walker to pull yourself up from a seated position. The knee walker is intended for use with a person in place on the walker. Use caution when moving from carpeted to hard surface floors and on any

Carex knee walkers offer mobility for those suffering from lower leg injuries and ailments who want to keep an active lifestyle.

- Knocks down quickly for easy transport and storage.
- Height adjustable handle bars and seat.
- Sturdy frame and tight turning radius allow for better maneuverability.

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