


walker instructions



10211-4 Shown
For use with all
Two Button, Deluxe Folding Walkers.



WEIGHT CAPACITY: 350 LBS.

drive™
www.drivemedical.com

Rev. 3.04.3015

wheelchair to walker transfer



Making certain wheels on the wheelchair are locked lift your body by pushing down on the armrests and carefully transfer one hand at time to the walker.



wheelchair to walker transfer




(with lower extremity injury to right leg or foot)

Making certain wheels on the wheelchair are locked, lift your body by pushing down on the armrests and carefully transferring first the left hand, then the right hand to the walker. (If injury to left leg or foot transfer the right hand first, followed by left.)



three-point gait

(for walking with walker)



Step 1

Align the middle of your foot with the back legs of the walker.

Step 2

Slide the walker forward and place it comfortably in front of you.

Step 3

Carefully, walk forward one step at a time.

0 0

custom walker adjustments


With your arm extended down, handle should be at wrist height.

With your hand placed on the hand-grip, your elbow should be flexed 20°-30°.

CAUTION:

Be aware of the following:

- Electrical cords in way
- Liquids on floor
- Throw rugs
- Toys or debris
- Furniture
- Doors and door steps
- Gaps in floor



WARNING:

1. Use wheels on front legs only or on all four legs with appropriate rear wheel brakes.
2. Do not push one side of the walker ahead and then other side without lifting the walker first.
3. Use Drive accessories only. Use of other manufacturers' accessories can result in improper fit and risk fall and/or injury.
4. Do not use on stairs or escalators.

lifetime limited warranty

Your Drive branded product is warranted to be free of defects in materials and workmanship for the lifetime of the product for the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Lifetime Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.

drive™

© 2015 Medical Depot, Inc. All rights reserved.
Drive is a trademark of Medical Depot, Inc.
Port Washington N.Y. 11050 USA Made in China

deluxe folding walker,
two button
with 5" wheels

item # 10211-4



EU Authorized Representative
Drive Medical Ltd, Fabbro
Elford, West Yorkshire,
United Kingdom WF5 9DP

www.drivemedical.com

drive™