

This information is provided for our Customers to ensure that they understand more about the products we have provided - both how they benefit the patient and how they can be safely used. If you have any questions about why this product was recommended for you or questions about a therapy regimen, please contact your physician or therapist.

If you have any questions or concerns about the product provided to you, please contact our office.

ALTERNATING PRESSURE PADS

INTRODUCTION to ALTERNATING PRESSURE PADS

Alternating Pressure Pads are designed to automatically change the pressure points beneath a patient every two and a half to four minutes. This is accomplished by alternately filling and emptying adjacent cells in the pad. Constantly changing pressure greatly reduces the danger of decubitus ulcers or pressure sores by varying the amount of contact pressure against the skin. This allows for better blood flow and improved tissue health.

MAINTENANCE & INSPECTION

The caregiver should check for the fill and empty movement daily by placing his or her open hand on the pad for a period of time long enough to allow the system to cycle from full to empty.

The family or caregiver should also inspect the tubing from the pad to the pump daily. It should not be twisted, pinched or kinked. Particular attention should be given to the ends of the tubes where they attach to the pad. This area is usually hidden by the sheet and is the most likely location of restricted air flow.

An overall visual inspection of the pad should be made each time the bed linens are changed. If the pad is ever removed from the bed for any reason, care should be exercised to place it back on the bed with the proper side up. Pads with a "top" and "bottom" side will be marked "This Side Up."

- No pins of any kind should be used in the bedding.
- Hot objects such as heating pads should not be placed on the alternating pressure pad.
- Heated under-blankets should not be used with the alternating pressure pad system.
- Solvents such as alcohol should not come into contact with the surface of the pad. A mild detergent and water should be used for cleaning.
- Caution must be exercised when providing bed care such as bathing the patient, so as not to expose the pump unit to accidental spills. The pump unit should not be exposed to liquids or moisture from open windows, aerosols or any other source.
- The power cord should not be allowed to come into contact with hot surfaces such as a heat register, and it should always be kept clear of the moving parts of the hospital bed.
- If unusual noises are heard coming from the pump unit, or if the pressure is not alternating properly, please discontinue use and call our office immediately.

Customers who qualify for a Support Surface may also need and qualify for a skin protection product or a combination cushion. Many also may require a Wheelchair, a Hospital-style Bed and mattress, Oxygen, Negative Pressure Wound therapy, or other supplies and products which help with the activities of daily living.

These may include:	Bedside Commode	Ostomy or Urological supplies
	Enteral Nutrition Supplies	Diabetic Supplies
	Protective Bed Pads	Patient Lifts & Ramp Systems

For more information about access to these products, please talk to our Representative.