

Nisus Negative Pressure Wound Therapy Pump

The Cork Nisus Negative Pressure Wound Therapy Pump is designed to provide gentle powered suction to treat and promote wound healing by aiding in the removal of excess exudates, infectious material, and tissue debris.



Latch design provides streamlined, one-handed canister release and removal



User friendly interface with easy to read LCD screen



FEATURES OF THE NISUS NEGATIVE PRESSURE WOUND THERAPY PUMP

The Cork Medical Products Nisus Negative Pressure Wound Therapy Pump is indicated for use in patients who would benefit from negative pressure wound therapy particularly as the device may promote wound healing by the removal of excess exudates, infectious material, and tissue debris.

- Recommended for patients with chronic, acute, traumatic, sub-acute and dehisced wounds, partial-thickness burns, pressure ulcers, diabetic ulcers, and venous ulcers, flaps, and grafts
- FDA approved for use in hospitals, long term care facilities, and home-care settings
- Light-weight and portable, weighing less than 2 lbs. with up to 48 hours of battery life
- Sleep feature technology allows user to dim light for longer lasting battery life and discreetness
- Step-by-step troubleshooting guide for alarms such as Canister Full, Vacuum Leak, and Blockage
- Multiple features for clinical enhancement such as: patient runtime, alarm log, lock restriction, Variable Intermittent and Continuous Therapy Modes

DEVICE SPECS

The Nisus Negative Pressure Wound Therapy Pump is part of a complete wound care system. Two essential components that actively work together to promote wound healing are the Nisus NPWT Canisters and the Cork NPWT Wound Dressing Kit. Please use under the advisement of a physician and carefully read all therapy and device instructions and safety information prior to use. Use only Cork Medical Products components with the Nisus Negative Pressure Wound Therapy Pump. Components are as follows:

| <u>Component</u> | <u>Model Number</u> |
|-------------------------|---|
| Nisus NPWT Pump | CMPP-100 |
| Nisus NPWT Canister | CPC-250 (250 mL) or CPC-500 (500 mL) |
| NPWT Wound Dressing Kit | WDK-100 (medium size) or WDK-200 (large size) |

- Battery charger
- Neoprene bag for carrying
- Medicare HCPCS Code: E2402