

# Patient Interfaces

## Roscoe Sapphire Series CPAP Masks

Roscoe's Sapphire CPAP Masks have an advanced mask design with a unique 3D cushion that forms an excellent seal and provides patients with superior comfort while protecting the skin from pressure points and irritations. It is both removable and replaceable. The innovative angle-adjustment forehead support has five levels of adjustment to meet a patient's individual needs and provide the best fit possible. The soft silicone pad contours to the forehead keeping the mask in place and stable throughout the night. The lightweight, quick-release headgear eliminates the need for constant adjustment. **Available in both full-face and nasal styles.**


Soft silicone pad conforms to forehead for increased stability.

Press the two clips together to adjust the angle of the forehead pad for an optimal fit.

Soft 3D cushion seal gently conforms to face and allows for jaw movement.

Quick-release headgear buckle

Built-in air vent disperses air quietly.



### Product Features:

- Lightweight advanced design ensures maximum comfort and compliance.
- Smaller facial footprint
- Comes complete with frame, 3D cushion, adjustable forehead pad, headgear, clips, and user manual with care instructions.
- Soft 3D cushion mask technology provides for an effective, yet comfortable seal. Cushions are removable and replacements are available.
- Innovative angle adjustable forehead pad has five levels of adjustment for the best fit possible.
- Quick-release headgear with Velcro® fasteners allows for easy removal and eliminates the need for constant adjustment.
- Quiet air vent minimizes noises so patients are not disturbed during the night.
- Available in full-face or nasal style in small, medium and large sizes
- Replacement seals available
- Latex-free
- 90 day limited warranty



Front



Back

**Attractive Packaging**  
**Easy-to-follow fitting instructions!**  
 Mask and headgear come complete in an attractive, resealable bag with easy fitting instructions on the back!





Sapphire Full  
Face Mask



Sapphire  
Nasal Mask

## Sapphire Series CPAP Masks

| Item #  | Description                                   | Qty  |
|---------|---|------|
| CPM-SNS | Sapphire Nasal Mask with Headgear, Small      | each |
| CPM-SNM | Sapphire Nasal Mask with Headgear, Medium     | each |
| CPM-SNL | Sapphire Nasal Mask with Headgear, Large      | each |
| CPM-SFS | Sapphire Full Face Mask with Headgear, Small  | each |
| CPM-SFM | Sapphire Full Face Mask with Headgear, Medium | each |
| CPM-SFL | Sapphire Full Face Mask with Headgear, Large  | each |

## Sapphire Series Replacement Parts

| Item # | Description                                 | Qty  |
|--------|---|------|
| 90605  | Sapphire Nasal Seal Replacement, Small      | each |
| 90607  | Sapphire Nasal Seal Replacement, Medium     | each |
| 90609  | Sapphire Nasal Seal Replacement, Large      | each |
| 90611  | Sapphire Full Face Replacement Seal, Small  | each |
| 90613  | Sapphire Full Face Replacement Seal, Medium | each |
| 90615  | Sapphire Full Face Replacement Seal, Large  | each |
| 90617  | Sapphire Series Replacement Headgear        | each |
| 90619  | Sapphire Series Replacement Headgear Clips  | pair |



*To Order:*  
 (800) 3-ROSCOE  
 fax: (440) 572-4261  
 roscoemedical.com